

Metabolomics Course

- Metabolomics --> why?
 - Data acquisition
 - How is the data?
 - Explain data parts (terms), the ones useful and add the visualization in the hands-on
 - Pre-processing steps: peak centroiding, peak detection, RT alignment, feature matching (steps in OpenMS and in xcms)
-
- **Add summaries from time to time**


Introduction to metabolomics

Complete set of small molecules (usually <1500 Da)
present in a biological system at a given time




Metabolome


Introduction to metabolomics




Introduction to metabolomics



Introduction to metabolomics



Introduction to metabolomics



Introduction to metabolomics

Applications


- 1 Discovery of new disease biomarkers
- 2 The understanding of molecular mechanisms in pathophysiological processes
- 3 Nutrition research
- 4 Drug discovery
- 5 Fermentation
- 6 Environmental Monitoring and Ecotoxicology

Introduction to metabolomics


Study types

- 1 Untargeted
- 2 Targeted
- 3 Semi-targeted

Metabolomics pipeline



Metabolomics pipeline



Metabolomics pipeline

Experimental design


Deal with confounding factors that could correlate or mask the biology under-study:

- Analytical: batch effects, injection order, etc.
- Biological: sex, age, body mass index, etc.


Biological samples randomization based on their physiological characteristics

Quality control (QC) samples are injected at the start, during and end of the acquisition to stabilize the signal and correct for technical variability.


Metabolomics pipeline




Metabolomics pipeline




Metabolomics pipeline



Metabolomics pipeline



Metabolomics data acquisition



Nuclear Magnetic Resonance (NMR)



Gas chromatography-Mass spectrometry (GC-MS)



Liquid chromatography-Mass spectrometry (LC-MS)

Metabolomics data acquisition

Nuclear Magnetic Resonance (NMR)

STRENGTHS

- 1 Non-destructive
- 2 Easily quantifiable
- 3 No prior separation required
- 4 Identification of novel compounds
- 5 High reproducibility

LIMITATIONS

- 1 Low-sensitivity
- 2 High detection limit
- 3 Large sample volumes

Metabolomics data acquisition

Gas chromatography-Mass spectrometry (GC-MS)

STRENGTHS

- 1 More mature and cost-effective
- 2 Excellent separation reproducibility
- 3 Availability of spectral libraries facilitates metabolite identification

LIMITATIONS

- 1 Sample derivatization required
- 2 Smaller coverage
- 3 Longer experimental times

Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)

STRENGTHS

- 1 Broader range of metabolites
- 2 Simpler sample preparation
- 3 High sensitivity


LIMITATIONS

- 1 Higher cost
- 2 Lower reproducibility
- 3 Reduced compatibility with volatile compounds

Preferred choice in metabolomics: **more than 70%** of the metabolomics studies were using LC-MS in 2022 .


Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)





Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)





Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)




Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)






Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)







Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)







Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)



Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)



Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)

Data Dependent Acquisition (DDA)

Ions are selected based on intensity.

The MS makes real-time decisions during the run: MS¹ scan → Selection of Top K most intense ions

Pros

High-quality MS/MS spectra

Cons

Biased toward abundant metabolites

Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)

Data Dependent Acquisition (DDA)

Ions are selected based on intensity.

The MS makes real-time decisions during the run: MS¹ scan → Selection of Top K most intense ions

Pros

High-quality MS/MS spectra

Cons

Biased toward abundant metabolites

Data Independent Acquisition (DIA)

It fragments all ions within predefined m/z windows.

Pros

Much more comprehensive coverage

Cons

MS/MS spectra are complex and overlapping

Metabolomics data acquisition

Liquid chromatography-Mass spectrometry (LC-MS)

Data Dependent Acquisition (DDA)

Ions are selected based on intensity.

The MS makes real-time decisions during the run: MS¹ scan → Selection of Top K most intense ions

Pros

High-quality MS/MS spectra

Cons

Biased toward abundant metabolites

Data Independent Acquisition (DIA)

It fragments all ions within predefined m/z windows.


Pros

Much more comprehensive coverage


Cons

MS/MS spectra are complex and overlapping


LC-MS data



LC-MS data




Computational metabolomics pipeline



Pre-processing

The goal is to convert **raw LC-MS data** into a **table of n peaks** defined by m/z , rt and an intensity value for each sample k .




— Pre-processing →

m/z	rt	I_1	...	I_k
mz_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...
mz_n	rt_n	I_{n1}	...	I_{nk}


Pre-processing

Metaboigniter pipeline


Pre-Processing



Identification (optional)




Predict metabolite formulas and structures with SIRIUS and CSI:FingerID.



Detect analogue molecules from your MS2 spectra with the machine learning tool MS2Query.

Feature Matrix



Features Intensities
Meta Data Identifications
Samples

Feature Linking

pyOpenMS

OPENMS

Pre-processing

Centroiding


- The MS has finite resolving power → ions are detected as profile peaks around their true m/z rather than as a single exact m/z
- First step: centroiding estimates the most probable peak center (Continuous mass spectra → series of discrete points)
- Reduction of amount of data without much loss of information.
- The next processing steps require centroided data.
- Many manufacturers already apply centroiding of the profile data (during acquisition or immediately after).

Pre-processing

Centroiding with OpenMS PeakPickerHiRes


Profile mode

MS1 spectrum (scan at RT = 5.2 min)



Centroid mode

MS1 spectrum (scan at RT = 5.2 min)




An Automated Pipeline for High-Throughput Label-Free Quantitative Proteomics

Pre-processing

Centroiding with OpenMS PeakPickerHiRes

Profile mode

MS1 spectrum (scan at RT = 5.2 min)




1. Detect regions of neighbouring data points that could be a peak. At least 3 points that satisfy:
 - Signal-to-noise filter (user-defined)
 - Evenly spaced across m/z
 - Peak shape: intensities first increase and then decrease

Pre-processing

Centroiding with OpenMS PeakPickerHiRes

Profile mode

MS1 spectrum (scan at RT = 5.2 min)




1. Detect regions of neighbouring data points that could be a peak. At least 3 points that satisfy:
 - Signal-to-noise filter (user-defined)
 - Evenly spaced across m/z
 - Peak shape: intensities first increase and then decrease

Pre-processing

Centroiding with OpenMS PeakPickerHiRes

Profile mode

MS1 spectrum (scan at RT = 5.2 min)




2. A cubic spline is fitted through the peak points to obtain a smooth continuous representation of the signal.

Pre-processing

Centroiding with OpenMS PeakPickerHiRes

Profile mode

MS1 spectrum (scan at RT = 5.2 min)




2. A cubic spline is fitted through the peak points to obtain a smooth continuous representation of the signal.

Pre-processing

Centroiding with OpenMS PeakPickerHiRes

Profile mode

MS1 spectrum (scan at RT = 5.2 min)



3. Find the peak region's apex (centroid):

- Derivative of the spline: $f'(x) = I'(x) = 0$
- Bisection method
- Peak region:

$$m/z = [100.001, 100.005]$$

$$f(100.001) = 4 \text{ and } f(100.005) = -2$$

$$m = (100.001 + 100.005) / 2 = 100.003$$

$$f(100.003) = +1 \text{ -- Apex in } [100.003, 100.005]$$




Centroid m/z -- $I(m/z)$ = centroid intensity

Pre-processing

Peak detection


Profile data

20171016_POOL_POS_3_105-134.mzML



Centroid data


20171016_POOL_POS_3_105-134.mzML



The **goal** is to identify the **chromatographic peaks** (peaks across the retention time axis) and return a **peak feature** defined by m/z, rt and intensity per sample.

Pre-processing

Peak detection with OpenMS FeatureFinderMetabo



1. Mass trace detection
2. Feature assembly

Kenar, E et al. (2014). Automated label-free quantification of metabolites from liquid chromatography-mass spectrometry data. Molecular & cellular proteomics.

Pre-processing

Peak detection with OpenMS FeatureFinderMetabo

1. Mass trace detection

$P = \{p_k\}$ where $p_k = \{rt_k, m/z_k, I_k\}$

- Sort P by decreasing intensity
- Remove peaks with $I_k <$ user-defined threshold



Filtered list P'

Mass trace $T \rightarrow$ List of n peaks $p_k \in P'$ with similar m/z and occur in adjacent scans of an LC-MS run.

Kenar, E et al. (2014). Automated label-free quantification of metabolites from liquid chromatography-mass spectrometry data. Molecular & cellular proteomics.

Pre-processing – Peak detection

OpenMS FeatureFinderMetabo

1. Mass trace detection

Mass trace T $\rightarrow$ List of n peaks $p_k \in P'$ with similar m/z and occur in adjacent scans of an LC-MS run.

1. A peak p_k is considered the start of a mass trace T : $T = \{p_k\}$
2. Extend trace T in time (both forward and backward):
 - Peaks in neighbouring scans
 - With Similar m/z values (probabilistic)
 - Add peaks to trace T
3. Trace extension stops when no more peaks can be added
4. Used peaks are removed to avoid duplication

Kenar, E et al. (2014). Automated label-free quantification of metabolites from liquid chromatography-mass spectrometry data. Molecular & cellular proteomics.


Pre-processing – Peak detection

OpenMS FeatureFinderMetabo

1. Mass trace detection

Two compounds can have:

- Almost identical m/z
 - Overlap in RT
-



1. Smooth the signal with LOWESS
2. Find:
 - Local maxima (peaks)
 - Minima between maxima (separation points)
3. Each split becomes a final mass trace

Kenar, E et al. (2014). Automated label-free quantification of metabolites from liquid chromatography-mass spectrometry data. Molecular & cellular proteomics.

Pre-processing – Peak detection


OpenMS FeatureFinderMetabo

2. Feature assembly

The goal is to find sets of isotopic mass traces that might originate from the same analyte.

A set is considered a valid isotope pattern if:

- Mass traces co-elute
- Correct m/z distances with respect to charge z
- Correct isotopic abundance ratios



Kenar, E et al. (2014). Automated label-free quantification of metabolites from liquid chromatography-mass spectrometry data. Molecular & cellular proteomics.

Pre-processing – Peak alignment and matching

The goal is to match features across different samples together:

1. Peak alignment across the retention time axis
2. Grouping corresponding features in multiple runs

Pre-processing – Hands on

Click on: Open in GitHub Codespaces

Mass spectrometry-based
Metabolomics introduction
documentation



Search

% + K

Metaboigniter

MS-based metabolomics

Metaboigniter



Metaboigniter

Open GitHub codespace

Use the following link to open a GitHub codespace with most of the required software installed:

⚠ If you do it manually, make sure to select the bigger machine with 4 cores and 16GB RAM



Open in GitHub Codespaces

Alternative setup of codespace

1

2


The screenshot shows the GitHub repository page for `biosustain/dsp_course_metabolomics_intro`. The repository is public and has 1 fork and 0 stars. A blue banner at the top indicates that the main branch is not protected. The repository contains several files and folders, including `.devcontainer`, `.github/workflows`, `.vscode`, `bin`, `data`, `material`, `results_prepared`, and `.gitignore`. The 'Code' button is highlighted with a green arrow, and the 'Codespaces' dropdown menu is open, showing options like 'New with options...', 'Configure dev container', 'Set up prebuilds', 'Manage codespaces', 'Share a deep link', and 'What are codespaces?'. A second green arrow points to the 'Codespaces' dropdown menu.

Codespace with 4 cores (and 16GB of memory)

- Not Safari
- Ignore Errors (todo tree)

Create codespace for
biosustain/dsp_course_metabolomics_intro

Branch This branch will be checked out on creation	 main ▾
Dev container configuration Your codespace will use this configuration	nextflow-training ▾
Region Your codespace will run in the selected region	Europe West ▾
Machine type Resources for your codespace	4-core ▾
Create codespace	



 **Your main branch isn't protected** Dismiss Protect this branch
Protect this branch from force pushing or deletion, or require status checks before merging. [View documentation](#).

 **main**  **4 Branches**  **0 Tags** T Add file  <> Code

 enryH and feliciaschulz Data analysis (#8) 6cd94bb · 48 minutes ago 9 Commits
 .devcontainer  set latest dockerfile (#4) 4 days ago
 .github/workflows  add docker image similar to proteomics course (#2) last week
 .vscode  set latest dockerfile (#4) 4 days ago
 bin Data analysis (#8) 48 minutes ago
 data Data analysis (#8) 48 minutes ago
 material Fix docs (#7) 11 hours ago

About

About Introduction to massspectrometry-based metabolomics (raw data processing and downstream analysis).

biosustain.github.io/dsp_course_...

-  Readme
-  Activity
-  Custom properties
-  1 star
-  0 watching
-  1 fork
- Report repository

Releases

Go to website

Repo Overview

Feel free to ask questions at
any point

enryH and feliciaschulz Data analysis (#8) 6cd94bb · 49 minutes ago 🕒 9 Commits		
📁 .devcontainer	🔧 set latest dockerfile (#4)	4 days ago
📁 .github/workflows	🔧 add docker image similar to proteomics course (#2)	last week
📁 .vscode	🔧 set latest dockerfile (#4)	4 days ago
📁 bin	Data analysis (#8)	49 minutes ago
📁 data	Data analysis (#8)	49 minutes ago
📁 material	Fix docs (#7)	11 hours ago
📁 results_prepared	Data analysis (#8)	49 minutes ago
📄 .gitignore	🔧 add docker image similar to proteomics course (#2)	last week
📄 2_data_analysis.ipynb	Data analysis (#8)	49 minutes ago
📄 2_data_analysis.py	Data analysis (#8)	49 minutes ago
📄 README.md	Updated Agenda and Location June 2026 (#6)	4 days ago
📄 conf.py	Data analysis (#8)	49 minutes ago
📄 index.md	Data analysis (#8)	49 minutes ago
📄 nextflow.config	🔧 set latest dockerfile (#4)	4 days ago
📄 requirements.txt	Data analysis (#8)	49 minutes ago
📄 requirements_docs.txt	🚀 initialize relatively bare bone repo for metabolomics cou...	3 months ago

Nextflow and nf-core pipeline template



Nextflow is a workflow executor
(analysis steps combined in an
execution graph)

- Reproducible and scalable

Nf-core is a collection of
workflows maintained by
nextflow and an open-source
community



A global community effort to collect a curated set of
open-source analysis pipelines built using Nextflow.

Pipelines

Browse the 130 pipelines that are currently available as part of nf-core.

Released 80

Under development 38

Archived 12

⌵ Last release ▾



proteinfamilies ✓

☆ 14

Generation and update of protein families

metagenomics protein-families proteomics

📦 1.1.0 released 4 days ago

epitopeprediction ✓

☆ 44

A bioinformatics best-practice analysis pipeline for epitope prediction
and annotation

epitope epitope-prediction mhc-binding-prediction


📦 3.0.0 released 5 days ago



Works with

- Linux machines (on servers)
- on MacOS



INPUTS

Metaboigniter – Metro Map





(...)



Samplesheet

sample	level	type	msfile
MS_A_POS	MS1	normal	data/MTBLS8735/MS_A_POS.mzML
MS_B_POS	MS1	normal	data/MTBLS8735/MS_B_POS.mzML
MS_C_POS	MS1	normal	data/MTBLS8735/MS_C_POS.mzML
MS_D_POS	MS1	normal	data/MTBLS8735/MS_D_POS.mzML
MS_E_POS	MS1	normal	data/MTBLS8735/MS_E_POS.mzML
MS_F_POS	MS1	normal	data/MTBLS8735/MS_F_POS.mzML
MS_QC_POOL_1_POS	MS1	normal	data/MTBLS8735/MS_QC_POOL_1_POS.mzML
MS_QC_POOL_2_POS	MS1	normal	data/MTBLS8735/MS_QC_POOL_2_POS.mzML
MS_QC_POOL_3_POS	MS1	normal	data/MTBLS8735/MS_QC_POOL_3_POS.mzML
MS_QC_POOL_4_POS	MS1	normal	data/MTBLS8735/MS_QC_POOL_4_POS.mzML
MSMS_2_A_CE20_POS	MS2	normal	data/MTBLS8735/MSMS_2_A_CE20_POS.mzML
MSMS_2_A_CE30_POS	MS2	normal	data/MTBLS8735/MSMS_2_A_CE30_POS.mzML
MSMS_2_A_CES_POS	MS2	normal	data/MTBLS8735/MSMS_2_A_CES_POS.mzML
MSMS_2_E_CE20_POS	MS2	normal	data/MTBLS8735/MSMS_2_E_CE20_POS.mzML
MSMS_2_E_CE30_POS	MS2	normal	data/MTBLS8735/MSMS_2_E_CE30_POS.mzML
MSMS_2_E_CES_POS	MS2	normal	data/MTBLS8735/MSMS_2_E_CES_POS.mzML

Spectra files


- Open standard, text based:
 - mzML files (indexed)
- Needs to be created from .d (bruker), .raw (Thermo), etc. vendor specific file formats
 - Conversion not (yet) integrated into workflow

One MS1 spectrum (DDA)


```
<spectrum id="controllerType=0 controllerNumber=1 scan=6" index="5" defaultArrayLength="736">
  <cvParam cvRef="MS" accession="MS:1000511" value="1" name="ms level" />
  <cvParam cvRef="MS" accession="MS:1000579" value="" name="MS1 spectrum" />
  <cvParam cvRef="MS" accession="MS:1000130" value="" name="positive scan" />
  <cvParam cvRef="MS" accession="MS:1000285" value="15248891" name="total ion current" />
  <cvParam cvRef="MS" accession="MS:1000127" value="" name="centroid spectrum" />
  <cvParam cvRef="MS" accession="MS:1000504" value="401.923461914063" name="base peak m/z" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
  <cvParam cvRef="MS" accession="MS:1000505" value="2704116.25" name="base peak intensity" unitAccession="MS:1000131" unitName="number of detector counts" unitCvRef="MS" />
  <cvParam cvRef="MS" accession="MS:1000528" value="375.884124755859" name="lowest observed m/z" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
  <cvParam cvRef="MS" accession="MS:1000527" value="1665.40612792969" name="highest observed m/z" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
  <scanList count="1">
    <cvParam cvRef="MS" accession="MS:1000795" value="" name="no combination" />
    <scan instrumentConfigurationRef="IC1">
      <cvParam cvRef="MS" accession="MS:1000016" value="6.0323342" name="scan start time" unitAccession="U0:0000031" unitName="minute" unitCvRef="U0" />
      <cvParam cvRef="MS" accession="MS:1000512" value="FTMS + p NSI Full ms [375.0000-1800.0000]" name="filter string" />
      <cvParam cvRef="MS" accession="MS:1000927" value="50" name="ion injection time" unitAccession="U0:0000028" unitName="millisecond" unitCvRef="U0" />
      <scanWindowList count="1">
        <scanWindow>
          <cvParam cvRef="MS" accession="MS:1000501" value="375" name="scan window lower limit" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
          <cvParam cvRef="MS" accession="MS:1000500" value="1800" name="scan window upper limit" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
        </scanWindow>
      </scanWindowList>
    </scan>
  </scanList>
  <binaryDataArrayList count="2">
    <binaryDataArray encodedLength="3152">
      <cvParam cvRef="MS" accession="MS:1000514" value="" name="m/z array" unitAccession="MS:1000040" unitName="m/z" unitCvRef="MS" />
      <cvParam cvRef="MS" accession="MS:1000523" value="" name="64-bit float" />
      <cvParam cvRef="MS" accession="MS:1000574" value="" name="zlib compression" />
      <binary>eJwT03LYXUawPHjVumjGA1lG5dTGtMi3MYUG32496hJ0iNX3EvLcjlNsEDNckGWIyoGxaIwiU3AsUwc22TJEnEm5Wimg48C466j15M12ZRgU1paCTPP+/399Xne9fceFk3TzIeylxna//0Vjck5ovM26h1oDLdI
    </binaryDataArray>
    <binaryDataArray encodedLength="4356">
      <cvParam cvRef="MS" accession="MS:1000515" value="" name="intensity array" unitAccession="MS:1000131" unitName="number of counts" unitCvRef="MS" />
      <cvParam cvRef="MS" accession="MS:1000523" value="" name="64-bit float" />
      <cvParam cvRef="MS" accession="MS:1000574" value="" name="zlib compression" />
      <binary>eJwTWHlCt+kXvILFSI2MTJi5ZWiaTTGdSl+EkT3KVtypUVl+dhR5TG4qo7Jk8FPWm32pjC1kuxU1Eu0JuSVDpozK0pdKpj3P/PV+3vu+73nPec45zznvFQRBUrdyQxIEQWzvf6ppFBz63sIYceZq06jafJHbNGr;
    </binaryDataArray>
  </binaryDataArrayList>
</spectrum>
```

Metaboigniter - Overview


Pre-Processing




Identification (optional)



Predict metabolite formulas and structures with SIRIUS and CSI:FingerID.



Detect analogue molecules from your MS2 spectra with the machine learning tool MS2Query.



Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}

Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}

Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}

Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}

Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}

Data analysis - Filtering and normalization

Input

Peak	m/z	rt	I_1	...	I_k
P_1	m/z_1	rt_1	I_{11}	...	I_{1k}
...	...	...	...	...	...
P_n	m/z_n	rt_n	I_{n1}	...	I_{nk}



Missingness due to errors or real
absence of the molecule

Data analysis - Filtering and normalization

Filtering using the 80% rules

Peak	m/z	rt	missigness (%)	l_1	...	l_k
P_1	m/z_1	rt_1	5%	l_{11}	...	l_{1k}
P_2	m/z_2	rt_2	8%	l_{21}	...	l_{2k}
P_3	m/z_3	rt_3	22%	l_{31}	...	l_{3k}
...	...	...	...	...	...	...
P_n	m/z_n	rt_n	42%	l_{n1}	...	l_{nk}

Data analysis - Filtering and normalization

Filtering using the 80% rules

Peak	m/z	rt	missigness (%)	I_1	...	I_k
P_1	m/z_1	rt_1	5%	I_{11}	...	I_{1k}
P_2	m/z_2	rt_2	8%	I_{21}	...	I_{2k}
P_3	m/z_3	rt_3	22%	I_{31}	...	I_{3k}
...	...	...	...	...	...	...
P_n	m/z_n	rt_n	42%	I_{n1}	...	I_{nk}

Data analysis - Statistical Analysis

Statistical analysis

MTBLS8735 - subset

Serum from

- three cardiovascular disease (CVD)
- Three healthy controls (CTR)
- Four quality control (QC) samples – pooled serum

<https://rformassspectrometry.github.io/Metabonaut/articles/dataset-investigation.html>

Data analysis – Hands on

Metabolomics Interpretation

Peak annotation vs peak identification

Identification levels according to the Chemical Analysis Working Group of the Metabolomics Standards Initiative:

Level 1 – **Identified** compounds require comparison of at least two independent properties (e.g., retention time and MS/MS spectra) with those of an authentic standard analysed under identical conditions.

Metabolomics Interpretation

Peak annotation vs peak identification

Identification levels according to the Chemical Analysis Working Group of the **Metabolomics Standards Initiative**:

Level 1 – **Identified** compounds require comparison of at least two independent properties (e.g., retention time and MS/MS spectra) with those of an authentic standard analysed under identical conditions.

Level 2 – Putatively **annotated** compounds based on similarity to MS/MS data in commercial or public databases (reference standards are not available).

Metabolomics Interpretation

Peak annotation vs peak identification

Identification levels according to the Chemical Analysis Working Group of the Metabolomics Standards Initiative:

Level 1 – **Identified** compounds require comparison of at least two independent properties (e.g., retention time and MS/MS spectra) with those of an authentic standard analysed under identical conditions.

Level 2 – Putatively **annotated** compounds based on similarity to MS/MS data in commercial or public databases (reference standards are not available).

Level 3 – Putatively characterized **compound class**.

Level 4 – **Unknown** feature.

Metabolomics Interpretation

Peak annotation vs peak identification

Identification levels according to the Chemical Analysis Working Group of the Metabolomics Standards Initiative:

Level 1 – **Identified** compounds require comparison of at least two independent properties (e.g., retention time and MS/MS spectra) with those of an authentic standard analysed under identical conditions.

Level 2 – Putatively **annotated** compounds based on similarity to MS/MS data in commercial or public databases (reference standards are not available).

Level 3 – Putatively characterized **compound class**.

Level 4 – **Unknown** feature.

LC-MS¹


Metabolomics Interpretation

Peak annotation using LC-MS¹

Peak annotation is a **major bottleneck** and direct matching of masses is **not enough**:

1

A single metabolite produces multiple peaks due to adducts and fragments.



C07481

Caffeine

Exact mass = 194.0804

— LC-MS + —→

$[M+H]^+$

$[M+Na]^+$

$[M+K]^+$

$[M+NH_4]^+$

Metabolomics Interpretation

Peak annotation using LC-MS¹

Peak annotation is a **major bottleneck** and direct matching of masses is **not enough**:

2

DB matching results in hundreds of potential candidates.



Metabolomics Interpretation

Peak annotation using LC-MS¹

Peak annotation is a **major bottleneck** and direct matching of masses is **not enough**:

1

A single metabolite produces multiple peaks due to adducts and fragments.

2

DB matching results in hundreds of potential candidates.

3

The metabolome is not completely annotated.



Computational algorithms that use m/z , rt and intensity to return putative annotations.

Metabolomics Interpretation

Peak annotation using LC-MS¹

Peak annotation is a **major bottleneck** and direct matching of masses is **not enough**:

1

A single metabolite produces multiple peaks due to adducts and fragments.

2

DB matching results in hundreds of potential candidates.

3

The metabolome is not completely annotated.



Computational algorithms



Low confidence level

putative annotations.

Metabolomics Interpretation

Peak annotation vs peak identification

Identification levels according to the Chemical Analysis Working Group of the Metabolomics Standards Initiative:

Level 1 – **Identified** compounds require comparison of at least two independent properties (e.g., retention time and MS/MS spectra) with those of an authentic standard analysed under identical conditions.

Level 2 – Putatively **annotated** compounds based on similarity to MS/MS data in commercial or public databases (reference standards are not available).

LC-MS²

Level 3 – Putatively characterized **compound class**.

Level 4 – **Unknown** feature.

Metabolomics Interpretation

Peak annotation using LC-MS² spectra


1. Spectral library matching
2. Fragmentation tree methods
3. In-silico fragmentation

Add definitions, mention that in-silico spectra can then be used for spectral library matching, and that we focus on spectral library matching, maybe add here limitations of spectral library matching and then mention in-silico and analogue search

Metabolomics Interpretation

Peak annotation using LC-MS² spectra



Spectral library matching



Metabolomics Interpretation

Peak annotation using LC-MS² spectra


Spectral library matching




Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Spectral library matching




← Similarity →



Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Cosine similarity



Query spectrum


Reference spectrum

Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Cosine similarity

1. Find peaks that match: $|mz(p) - mz(p')| < t$ (e.g., 0.05 Da)




Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Cosine similarity

1. Find peaks that match: $|mz(p) - mz(p')| < t$ (e.g., 0.05 Da)



Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Cosine similarity

1. Find peaks that match: $|mz(p) - mz(p')| < t$ (e.g., 0.05 Da)
2. Intensities of matching peaks are converted to vectors

$I = [38, 46, 101]$

$I' = [30, 59, 85]$

Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Cosine similarity

1. Find peaks that match: $|mz(p) - mz(p')| < t$ (e.g., 0.05 Da)

$I = [38, 46, 101]$

2. Intensities of matching peaks are converted to vectors

3. Cosine similarity:

$I' = [30, 59, 85]$

$$\cos(\theta) = \frac{A \cdot B}{\|A\| \|B\|} = \frac{\sum_{i=1}^n A_i B_i}{\sqrt{\sum_{i=1}^n A_i^2} \sqrt{\sum_{i=1}^n B_i^2}} = 0.98$$

Metabolomics Interpretation

Peak annotation using LC-MS² spectra

What if ... Two molecules differ by a methyl group:

Molecule A precursor = 300 Da

Molecule B precursor = 314 Da



Difference = 14 Da

Many fragments may shift +14 Da

— Cosine similarity
approach —>

No matching peaks


Similarity = 0

But molecules are structurally very similar... Modified cosine similarity

Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Modified cosine similarity



Query spectrum

Reference spectrum


Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Modified cosine similarity

1. Find peaks that match: $|mz(p) + M - mz(p')| < t$


Precursor mass difference $M = PM(S') - PM(S) = 14$ Da



Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Modified cosine similarity



1. Find peaks that match: $|mz(p) + M - mz(p')| < t$
Precursor mass difference $M = PM(S') - PM(S) = 14$ Da
2. Intensities of matching peaks are converted to vectors
3. Cosine similarity:


$$\cos(\theta) = \frac{A \cdot B}{\|A\| \|B\|} = \frac{\sum_{i=1}^n A_i B_i}{\sqrt{\sum_{i=1}^n A_i^2} \sqrt{\sum_{i=1}^n B_i^2}}$$

Metabolomics Interpretation




Peak annotation using LC-MS² spectra

Spectral library matching

▲ Apratoxin A ($m/z = 840.497$)



◆ Apratoxin A - 30.010 Da ($m/z = 810.487$)



Metabolomics Interpretation

Peak annotation using LC-MS² spectra

Spectral library matching

However...


Small chemical modifications → **Low** modified cosine scores

Ideally, we want an approach that understands chemical structure!

Metabolomics Interpretation

Peak annotation using LC-MS² spectra


Machine learning-based similarities



Metabolomics Interpretation

Peak annotation using LC-MS² spectra


Machine learning-based similarities



Metabolomics Interpretation


Peak annotation using LC-MS² spectra

Machine learning-based similarities




Metabolomics Interpretation

Peak annotation with MS2Query




Metabolomics Interpretation

Pathway enrichment




Metabolomics Interpretation

Pathway enrichment




Metabolomics Interpretation

Pathway enrichment



Metabolomics Interpretation

Pathway enrichment



Metabolomics Interpretation


Over-Representation Analysis (ORA)

Test if the proportion of **differentially expressed molecules** within a particular pathway is statistically higher or lower than would be **expected by chance**.

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Test if the proportion of **differentially expressed molecules** within a particular pathway is statistically higher or lower than would be **expected by chance**.



Picart-Armada, Sergio, et al. "FELLA: an R package to enrich metabolomics data." BMC bioinformatics

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Background: 100 metabolites measured (that belong to 3 pathways)

Pathway	Total # of metabolites
Pathway A	10
Pathway B	20
Pathway C	30
Not in a pathway	40
Total	100

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Background: 100 metabolites measured (that belong to 3 pathways)

10 differentially expressed metabolites

Pathway	Total # of metabolites	Significant metabolites
Pathway A	10	5
Pathway B	20	2
Pathway C	30	1
Not in a pathway	40	2
Total	100	10

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Background: 100 metabolites measured (that belong to 3 pathways)

10 differentially expressed metabolites

How many metabolites would we expect by chance?

Pathway	Total # of metabolites	Significant metabolites	Random
Pathway A	10	5	$10 \times (10/100) = 1$
Pathway B	20	2	$10 \times (20/100) = 2$
Pathway C	30	1	$10 \times (30/100) = 3$
Not in a pathway	40	2	
Total	100	10	

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Background: 100 metabolites measured (that belong to 3 pathways)

10 differentially expressed metabolites

How many metabolites would we expect by chance?

Pathway	Total # of metabolites	Significant metabolites	Random
Pathway A	10	5	$10 \times (10/100) = 1$
Pathway B	20	2	$10 \times (20/100) = 2$
Pathway C	30	1	$10 \times (30/100) = 3$
Not in a pathway	40	2	
Total	100	10	

Metabolomics Interpretation

Over-Representation Analysis (ORA)

Contingency table for Pathway A

	In A	Not in A	Total
Significant	5	5	10
Not significant	5	85	90
Total	10	90	100

Hypergeometric distribution: $P(X \geq 5)$

Interpretation – Hands on

Acore enrichment functionality



Search % + K

Acore Documentation

API usage examples

- Metabolomics data filtering
- Normalization of samples
- Batch correction of samples
- Imputation data (MS-example)
- Imputation (Metabolomics example)
- Exploratory Analysis
- Metabolomics drift correction
- Differential regulation: T-test for two groups
- Differential regulation: ANOVA across more than two groups and posthoc t-tests
- Differential regulation: ANCOVA for two groups
- Enrichment analysis**



Enrichment analysis

Is done separately for up- and downregulated genes as it's assumed that biological processes are regulated in one direction.

► Show code cell source

	identifier	group1	group2	pvalue	padj	rejected	log2FC	FC
98	O43432	deceased	alive	0.000	0.001	True	-0.741	0.598
32	Q96JJ3	deceased	alive	0.001	0.018	True	-0.553	0.682
99	O43175	deceased	alive	0.000	0.001	True	1.321	2.498
97	P39059	deceased	alive	0.000	0.000	True	1.600	3.031

Running the enrichment analysis for the up- and down regulated protein groups separately with the default settings of the function, i.e. a log2 fold change cutoff of 1 and at least 2 protein groups detected in the set of proteins defining the functional annotation.

```
ret = acore.enrichment_analysis.run_up_down_regulation_enrichment(  
    regulation_data=diff_reg,  
    annotation=annotations,  
    pval_col="padj",  
    min_detected_in_set=2,  
    lfc_cutoff=1,  
)
```

No significant enrichment found with the given parameters. Returning an empty DataFrame.



Q Search

⌘ + K

ACore Documentation

API usage examples

Metabolomics data filtering

Normalization of samples

Batch correction of samples

Imputation data (MS-example)

Imputation (Metabolomics example)

Exploratory Analysis

Metabolomics drift correction

Differential regulation: T-test for two groups

[Differential regulation: ANOVA across more than two groups and posthoc t-tests](#)

Differential regulation: ANCOVA for two groups

Enrichment analysis

Reanalysis

[Data Analysis PXD040621](#)



► Show code cell source



Enrichment Analysis regulated vs non-regulated




Some slides from
transcriptomics course

Functional analysis for biological interpretation

- **Enrichment analysis**

- Knowledge-based approach for biological interpretation of omics data, e.g. RNAseq results
- **Goal:** find biological terms (curated gene sets for biological processes) that are enriched in the data for biological interpretation (no single-gene interpretation)
- Biological term/pathway lists are obtained from KEGG, GO etc or can be manually curated (often .gmt files)
- Different types of enrichment analysis



Functional analysis for biological interpretation

- **Gene-set enrichment analysis (GSEA)** – all genes ranked by fold-change and p-value

List of all genes

	Gene symbol	p-val	Differentially expressed	FC
1	<i>HAO1</i>	0.0001	Yes	-5.6
2	<i>RDH16</i>	0.00045	Yes	-7.6
3	<i>ALDOB</i>	0.00009	Yes	-6.7
4	<i>AKR1C4</i>	0.0051	Yes	-10.0
5	<i>ABCB11</i>	0.067	Yes	10.0
6	<i>BAAT</i>	0.0034	Yes	23.5
7	<i>SLC2A2</i>	0.00990	Yes	2
8	<i>FCER1G</i>	0.00673	Yes	6.7
...				
6532	<i>OSM</i>	0.00009	Yes	9
6533	<i>VEGFA</i>	0.067	Yes	8.5
6534	<i>ZAP70</i>	0.0034	Yes	7.8

$$\text{RANKING} = \text{SIGN}(\text{FC}) * -\log_{10}(\text{p-value})$$

Ranked list of genes

	Gene symbol	Ranking
1	<i>JAK9</i>	12.9
2	<i>RDH16</i>	11.4
3	<i>ALDOB</i>	11.3
5	<i>AKR1C4</i>	10.7
6	<i>ABCB11</i>	8.6
7	<i>BAAT</i>	3.6
8	<i>SLC2A2</i>	0.0099
9	<i>FCER1G</i>	0.00673
	...	
6532	<i>MyD88</i>	-.9
6533	<i>TRIF</i>	-9.37
6534	<i>IRAK1</i>	-15.6

SIGNIFICANT + UPREGULATED

NON-SIGNIFICANT

SIGNIFICANT + DOWNREGULATED

Functional analysis for biological interpretation

- **Gene-set enrichment analysis (GSEA)** – Which biological processes are enriched at the top and at the bottom of the ranked gene list?

	Gene symbol	Ranking
1	<i>PTPN22</i>	12.9
2	<i>RDH16</i>	11.4
3	<i>IL1B</i>	11.3
5	<i>TNF</i>	10.7
6	<i>ABCB11</i>	8.6
7	<i>SPON2</i>	3.6
8	<i>SIC2A2</i>	0.0099
9	<i>TRIM55</i>	0.00673
	...	
6532	<i>OSM</i>	-9
6533	<i>TLR6</i>	-9.37
6534	<i>ZAP70</i>	-15.6

RANKED LIST



(181 GENES)

interleukin-6 production
(GO:0032635)*

PTPN22, IL1B, CSK, TNF,
TRAF6, TLR2, RABGEF1,
PYCARD, SIRPA, SPON2,
TLR7, TLR6, TRIM55, ZG16,
TYROBP, UCN, XBP1...

GENE SET

*GO = Gene Ontology

Result will show genes and biological terms that are

Up-regulated

[illegible]

Down-regulated

[illegible]

Not enriched

[illegible]

Statistical metrics:

- Enrichment score (ES) of a biological term represents the running sum when walking down the list
- Normalised enrichment score (NES): ES is normalised to make biological terms with different number of genes comparable
- p-value from permutations (adjusted)

Functional analysis for biological interpretation

- **Enrichment analysis – details and summary**

	Over-representation analysis (ORA)	Gene set enrichment analysis (GSEA)
Input gene list	<u>Filtered list</u> , e.g. $p_{\text{adj}} < 0.05$ and $\log_2\text{FC} > 1 $ → Genes not passing the thresholds will be missed	<u>Entire gene list ranked</u> by fold-change (with sign) and/or p_{adj} → No pre-filtering of genes is an advantage
Separation of up- and down-regulated genes	Can be done	No
Core statistical result and tests used	• Typically an adjusted p value • Fisher's exact test (one-sided), hypergeometric distribution	• Normalized enrichment score (NES) and adjusted p value • Modified Kolmogorov-Smirnov and permutation test
Software packages (selection)	mulea, g:profiler, Enrichr, clusterProfiler, GSEApy (via Enrichr)	fgsea, clusterProfiler, GSEApy